EXP:	TEST	FOR	ANIO	NS
------	------	------------	------	----

NAME___

PERIOD

PROCEDURE: For each anion test, use a squirt of each chemical in a test tube. Record your observations after each addition of a chemical.

To TEST FOR	ADD	OBSERVATIONS		
GROUP 1 Precipitate with BaCl,				
1. Sulfate, SO ₄ -2	BaCl₂ and HCl	1.		
2. Carbonate, CO₃⁻²	BaCl₂ and HCl	2.		
3. Phosphate, PO ₄ -3 (If no reaction, go bac solution and do step 4	4.)	3 .		
GROUP 2 Precipitate	WILL AGIVO3			
4. Phosphate, PO ₄ -3	AgNO ₃	4.		
4a. PO ₄ -3, confirming tes	4a.			
5. Chloride, CI -1	AgNO ₃ , NH ₄ OH, HNO ₃	5.		
6. Bromide, Br ⁻¹	AgNO ₃ , NH ₄ OH, HNO ₃	6.		
6a. Br ⁻¹ , confirming test, starch solution, add 1	6a.			
7. lodide, I ⁻¹	AgNO ₃ , NH ₄ OH, HNO ₃	7 .		
7a. I ⁻¹ , confirming test, A starch solution, add 1	7a .			

GROUP 3 Soluble group, will not precipitate

8. Nitrate, NO₃-1

FeSO₄ and conc. H₂SO₄

8.

Use equal amounts of NaNO₃ and FeSO₄, then add a small amount of conc. H₂SO₄ to the tilted test tube, allowing

it to run down the side. Look for the brown ring.